DAILY CLEANING FOR TRANSFER TABLE VACUUM PLATE

Embrex® BioDevices vacuum plate cleaning relies on the Hatchery Sanitation Program to ensure proper elimination of microbial challenges, both bacterial and fungal.

NOTE: The probability of contamination increases when proper cleaning and storage steps ARE NOT followed.

PREPARATION AND OPERATION OF THE VACUUM HEAD CLEANING BATH:

1. At completion of the injection, after pressing an e-stop, disconnect the transfer table from BioDevice vaccinator:

- Disconnect the electrical cable(s) between the vaccinator cabinet and the transfer table.
- Disconnect the AC plug/cable between the vaccinator cabinet and the transfer table if a separate connector is used.
- Disconnect the pneumatic (air) connector between the vaccinator pneumatic cabinet and the transfer table.
- Release the mechanical clamps, allowing the vaccinator and transfer table frames to be separated.
- Move the transfer table away from the vaccinator before starting any vaccinator cleaning cycles.

2. Remove the vacuum plate from the transfer table.

NOTE: The transfer plate should be in the down position and resting on or just above the table surface for this step.

- Remove the internal/external stripper weights.
- Release the external clamps on the plate.
- Separate the “lower” vacuum plate (removable plate with vacuum cups) from the “upper” plate (plate attached to transfer cylinder).
- Raise the “upper” plate (by hand) and insert the hitch pins in the guide shafts to hold it in position.

3. Flush/rinse visible organic matter from both sides of the “lower” transfer plate

- Water from a hose outfitted with a simple garden hose sprayer is sufficient. A soft brush can be used in removing the organic matter, but care must be exercised to avoid damaging vacuum cups.

4. Place the vacuum transfer plate into the cleaning bath provided by Zoetis.

- Prepare the “CLEANING BATH SOLUTION” using 240 ml of a “NON-FOAMING” chlorinated cleaner approved by Zoetis.

NOTE: The chlorinated cleaner should be added after the cleaning tank is filled with water. It is highly recommended that warm water is utilized for cleaning, as efficiency is improved. Water should be between 100 F and 115 F.

- Since chemical makeup of chlorinated cleaners varies, consultation with a Zoetis field service representative is strongly recommended prior to purchase or use of a non-foaming chlorinated cleaner.

5. Fill cleaning bath with enough water to cover the transfer plate.

- Water level should be no more than 1 inch above the vacuum plate.

NOTE: The cleaning bath will perform best with warm water.

- “Freezing” cold water is not recommended and will have issues with efficacy of cleaning. Do not use water below 65 F for cleaning.
- Scalding hot water can damage the vacuum cups and be a safety hazard for the operator. Do not use water > 115 F.
- Check plate alignment to cleaning bath air manifold. Ensure that bubbles are rising through each vacuum cup fitting.
- Inspect and unclog the vacuum cup fittings as needed so that the bubbling action is able to clean inside of fittings.

NOTE: Allow cleaning bath to run for:

- 30 minutes when using warm water
- 1 hour when using cool water (see above....efficacy of cleaning is reduced and may not be adequate if water is too cold)
- Once allotted time has passed according to water temperature as mentioned above, drain water from cleaning bath and rinse/flush bath using tap water.
- Rinse transfer plate thoroughly with tap water after removal from cleaning bath.

NOTE: It is recommended that at least every two weeks the vacuum cups should be removed from the vacuum cup fittings and all vacuum cup fittings on the lower plate be individually inspected/cleaned with a mild detergent (1 ounce dishwashing detergent per gallon water) and small brush or similar to assure vacuum performance and removal of any organic debris in vacuum pathway.

6. Remove and clean the 3 transfer table vacuum hoses.

- Inspect hoses while they are disconnected and remove visible organic matter.
- Water from a hose outfitted with a simple garden hose sprayer is sufficient to remove the organic matter.
- Cover the exposed vacuum fittings with a small piece of vacuum hose to prevent water from entering the vacuum motor/mufflers.

7. Place the hoses in the cleaning bath.

- The hoses can be placed into the cleaning bath with the vacuum transfer head.
- Vacuum hoses should remain in the tank during cleaning vacuum transfer head (30 minutes for warm water or 1 hour for cool water).
- Remove the hoses after cleaning and flush with water, removing any residual chlorinated cleaner.
- Hang the hoses VERTICALLY to allow for drainage and drying.

NOTE: Water left in the hoses can enter the vacuum pump assembly and cause premature failure of the pump.

8. The “upper” plate (connected to transfer cylinder) should be cleaned using soap and water with a scrub brush or mild abrasive pad.

- Remove all organic matter from the transfer head “upper” plate surfaces.

9. The vacuum elbow, located on transfer head “upper” plate, should be cleaned using the 1-inch brush provided by Zoetis.

- The vacuum elbow can be removed to facilitate cleaning and closer inspection.
- Removal of all organic matter is the desired goal.

10. The vacuum selector switch is located below the transfer table and should be cleaned using the 1-inch brush provided by Zoetis

- Ensure that all organic matter is removed.

NOTE: Verify that the transfer table air is disconnected from the pneumatic cabinet prior to cleaning the vacuum selector switch.

In some instances, contact or exposure to organic matter and/or cleaning solutions may present health hazards. Appropriate safety equipment and procedures should be used at all times. This cleaning process may also generate wastes and byproducts that require proper handling and disposal.

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