DOCUMENT INFORMATION & PRINTING INSTRUCTIONS

ECN2981 – 6/19/2014 - JDC

NOTES:

- 1. THE FOLLOWING PAGES, NOT INCLUDING THIS COVER SHEET, TO BE PRINTED ON THE FRONT & BACK OF ONE PAGE, AND LAMINATED BASED ON THE FOLLOWING SPECIFICATIONS
- 2. ZOETIS PMS SPECIFIC 2-COLOR OR 4-COLOR PROCESS PRINTED ON 11” X 17” WITH SHORT-SIDE BLEEDS, 60# COVER-STOCK PAPER, DOUBLE-SIDED
- 3. LAMINATE INFORMATION: 5-MIL ENCAPSULATED WITH 1/4” BORDER, CORNERS ROUNDED
- 4. LABEL AND ARTWORK TO COMPLY WITH REQUIREMENTS OF ZOETIS LOGO STYLE GUIDE, 4-COLOR PMS MATCH OR PMS SPOT / GRADIENT COLORS

PROPRIETARY AND CONFIDENTIAL  
THE INFORMATION CONTAINED IN THIS DRAWING IS THE  
SOLE PROPERTY OF ZOETIS, INCORPORATED. ANY  
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN  
PERMISSION OF ZOETIS, INCORPORATED IS PROHIBITED.

KEY CONSIDERATIONS TO SUPPORT IN OVO VACCINE STERILITY

1. Vaccine Preparation and Handling

- **Maintain Bio-Security of the Vaccine preparation lab – restrict access, no through traffic, and no one allowed except vaccine prep personnel**
- **Consider installation of Laminar Flow hood for Inovo Vaccine preparation if Vaccine lab is compromised with regards to ability to disinfect/clean, air quality, shared work space (day of age spray vaccines).**
- **Strict asepsis during vaccine preparation.**
  - Wash hands with soap and water, dry thoroughly with clean paper towel prior to every mixing session. Hand sanitizer may be utilized after washing, but is not a replacement or alternate procedure for washing.
  - Use of disposable latex gloves is not an alternative for washing hands. Discourage use of gloves in mixing and handling processes as the gloves limit the ability to detect wet surfaces on vaccine ampules, vaccine bag, etc.
  - Spray and wipe down mixing area (lab counter top) with 70% isopropyl alcohol or dilute chlorine bleach prior to each mixing session (80 ml of 5% bleach per gallon water = 1000 ppm Cl).
  - Separate needle and syringe for each additive, single session usage, and correct sizes for use (18 g needle 1 and ½ inches long, 20 cc syringe). Utilize 20 cc syringe for diluent buffer (5-8 ml in syringe) pre-load prior to addition of vaccine from vial.
  - Quality control of diluent.
    - Store in clean, dry environment at room temperature (65-75 F).
    - Check clarity, color, small amount of air present, no distension of bag (bloating).
    - No condensation in outer wrap. Spray bag with disinfectant (alcohol or dilute chlorine) and dry thoroughly prior to use if condensation present.
  - Swab all needle entry ports (diluent, dye, lyophilized vaccine) with sterile alcohol swab before all needle entries.
  - **Dry ampules and hands thoroughly with new paper towel for each mixing episode prior to opening thawed (Marek’s vaccine) ampules.**
  - Never touch needle with finger.
  - Discard any questionable elements of vaccine prep, including vaccine if necessary, if mistakes are made.
  - **Common mistake is overflow of MD ampule during ampule rinsing procedure...discard rinse diluent, needle/syringe, (empty) ampules immediately if this occurs...DO NOT ADD DILUENT RINSE TO VACCINE BAG IF OVERFLOW OCCURS DURING RINSE PROCEDURE.**
  - Thaw bath should be clean and DRY during storage and prior to use. Clean with detergent, wipe/rinse clean with new paper towel, final spray with 70% isopropyl alcohol and allow to dry.
  - Include chlorine in thaw bath (15 ml of 5.0% bleach per gallon of distilled water results in ~ 200 ppm chlorine).
  - Do not use tap water in thaw bath, use distilled water (plus chlorine as per above). Add 15 ml of 5.0% chlorine to ice bath in transport cooler for chilling mixed vaccine prior to hanging on Inovoject system.
  - Clean and disinfect ice packs and MD ampule holder (alcohol or dilute chlorine spray) prior to use and at the end of each run day. Use only plastic or foil covered ice packs.
  - Clean and disinfect freezer and refrigerator minimum of monthly.
  - Maintain daily vaccine log – document vaccine/diluent serial numbers, time mixed, usage, etc.
  - Establish microbial sampling of mixed vaccine (every bag) on daily basis
- **Follow strict asepsis in attaching vaccine to Inovoject® system.**
  - Dry hands, vaccine bag and port entry area thoroughly with paper towel just prior to hanging vaccine on Inovoject system. Assure dry working area on vaccine bag at port entry point.
  - Do not remove cover of spike entry port (blue tag or twist port cover) on vaccine bag or saline (start up) until just prior to attaching to vaccine delivery system on machine.
  - Swab spike and entry port of vaccine bag with sterile alcohol swab prior to attachment to the vaccine delivery system.
- **Adhere strictly to daily, weekly, and monthly cleaning of vaccine preparation area (refer to Embrex Inovoject Process Cleaning Procedures).**
- **ONLY the vaccine mixing person is allowed to handle vaccine for in ovo use.**
- **Commit to a scaled microbial sampling program (daily, weekly, monthly regime) of vaccination process.**

2. Inovoject Process

- Adhere strictly to start up and shut down procedures for the Inovoject system as outlined in User’s Guide. Store Inovoject system in a dry, clean area overnight.
- Adhere to disposable change out schedule outlined by Zoetis.
- At shut down, be sure to clean all containers used on on-board cleaning cart (alcohol, sanitizer, cleaner, distilled water) daily and allow drying upside down overnight. Remember to add sufficient amount of alcohol to distilled water and alcohol container for overnight storage.
- Pay attention to general cleanliness of the Inovoject system, including drip pans, injection tooling, and vacuum cups on transfer head.
- Remove as many rotten, cracked, or broken eggs as possible prior to injection, collect and dispose appropriately.
- Monitor and confirm the mixing and proper use of Inovotabs™ as chlorine source for sanitation fluid.
- Monitor operation of QC plate, paying close attention to sanitation dispense (step 2 of process, after vaccine dispense verification).

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