DOCUMENT INFORMATION & PRINTING INSTRUCTIONS

ECN2981 – 6/19/2014 - JDC

NOTES:

1. THE FOLLOWING PAGES, NOT INCLUDING THIS COVER SHEET, TO BE PRINTED ON THE FRONT & BACK OF ONE PAGE, AND LAMINATED BASED ON THE FOLLOWING SPECIFICATIONS
2. ZOETIS PMS SPECIFIC 2-COLOR OR 4-COLOR PROCESS PRINTED ON 11” X 17” WITH SHORT-SIDE BLEEDS, 60# COVER-STOCK PAPER, DOUBLE-SIDED
3. LAMINATE INFORMATION: 5-MIL ENCAPSULATED WITH 1/4” BORDER, CORNERS ROUNDED
4. BRASS 1/2” GROMMETS WITH 1” BOOK RINGS TO ATTACH SPANISH/ENGLISH VERSIONS OF START-UP & SHUT-DOWN PROCESSES
5. LABEL AND ARTWORK TO COMPLY WITH REQUIREMENTS OF ZOETIS LOGO STYLE GUIDE, 4-COLOR PMS MATCH OR PMS SPOT / GRADIENT COLORS

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE
SOLE PROPERTY OF ZOETIS, INCORPORATED. ANY
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN
PERMISSION OF ZOETIS, INCORPORATED IS PROHIBITED.

EMBREX[®] INOVOJECT[®] BIODEVICES STARTUP PROCEDURE

EMBREX[®]

1. Machine Connections

- Attach vacuum pickup plate and stripper weights.
- Attach transfer table to the vaccinator and egg remover, if applicable.
- Mechanical
- Electrical
- Pneumatic

2. Prepare Chemicals

- Refer to user/operators manual or placard.

3. Start AM Cleaning Cycle

- Drop needles and inspect to ensure proper flow through **all** needles. ***NOTE – Pause cycle and correct flow if issues are observed.***
- When complete, the interface will read “AM Cleaning Cycle Complete.”
- Vaccine lines are now stored in alcohol and ready for saline and vaccine.

4. Sanitation System Setup

- Place sanitation pump tube in a clean container with 3 gallons (11.5 L) clean water.
- Press “**Sanitation Dispense**” to start sanitation pump and flush with 2 gallons of water.
- Place sanitation tube in prepared sanitizer.
- Press “**Sanitation Dispense**” to start sanitation pump and prime system with sanitizer for **at least** 2 minutes.
- Press “**Sanitation Dispense**” to stop sanitation pump after 2 minutes; the system is ready for operation.

5. Prime Embrex Inovoject with Sterile Water/Saline and Vaccine

- Drop needles.
- Clamp off the main yoke tube to vaccine manifold using forceps to avoid introducing air into the sterile system.
- Attach a bag of sterile saline/water to the vaccine system. **Swab spike with alcohol swab.**
- Prime vaccine delivery system with sterile saline/water.
- When finished, clamp off the vaccine line, remove prepared vaccine bag from transport cooler, dry prepared vaccine bag with fresh paper towel and connect to spike. **Swab spike with fresh alcohol swab before spiking vaccine and avoid allowing air into the system.**
- Prime vaccine delivery system until vaccine can be seen at every needle.
- With vaccine delivery system primed, place the vaccine yoke tubing in the vaccine sensor.
- Raise needles.
- Run the QC plate to verify complete priming of the vaccine delivery system. Verify all positions dispensing vaccine prior to egg injection.
- The BioDevice is ready to inject.

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EMBREX[®] INOVOJECT[®] BIODEVICES SHUTDOWN PROCEDURE

1. Complete QC plate check prior to starting the PM cleaning cycle to ensure the vaccine fluid pathways are clear.

2. Cleaning Pickup Plate

- Refer to placard for **Daily Transfer Table Vacuum Plate Cleaning**

3. PM Cleaning Cycle (refer to users/operators manual)

- E-stop the device and separate the transfer table from the vaccinator prior to starting the cycle.
- Connect spike to the cleaning cart line.
- Standard Inovoject – Vaccine Manifold fittings in “up” position and tech mate pump lid unlatched
- Start PM cleaning cycle
- Drop the needles to check for proper flow and bent needles. This should be completed in the first step of the cycle and the cycle should be paused to correct any issues that are observed.
- When complete, the interface will read “**PM Cleaning cycle complete.**”

WARNING: DO NOT ATTEMPT TO WASH ANY BioDevice WHILE ENERGIZED.

4. Sanitation System Shutdown

- Remove sanitation pump tube from the sanitizer container and properly discard remaining sanitizer.
- Place sanitation pump tube in a clean container with 3 gallons (11.5 L) clean water.
- Press “**Sanitation Dispense**” and flush with a minimum of 2 gallons (8.5 L) clean water.
- Once flush with water is complete, place sanitation tube in the alcohol storage container.
- Press “**Sanitation Dispense**” to start flush with alcohol for **2 minutes**. Press “**Sanitation Dispense**” to stop.
- The system is now stored in alcohol.
- Discard all remaining alcohol from alcohol storage container.
- Wash alcohol storage container with soap and water, rinse and dry.
- Add 0.5 gallons of fresh alcohol to alcohol storage container. The sanitation pump tube filter should be placed back into the alcohol container and remain submerged in alcohol for storage.

5. Clean the Transfer Table

- Rinse the machine with water. **No high-pressure washers.**
- Remove all vacuum hoses, rinse to remove debris and place in the cleaning tank with transfer head plate.
- Scrub all device surfaces with a scrubbing pad and antibacterial detergent (dishwashing liquid preferred).
- Key areas to scrub are the pickup plate, vacuum hose elbow and the pulleys beneath the table, then rinse

6. Clean the Vaccinator

- **After** the PM cleaning cycle has completed, **unplug the vaccinator.**
- Rinse machine with water. **No high-pressure washers.**
- **Use caution when rinsing. Avoid direct** spray of the interface (HMI), vaccine pump(s) and sensors.
- Scrub all surfaces with a scrubbing pad and antibacterial detergent (dish washing liquid).
- **Key areas to scrub are tooling and pulleys.**
- Rinse the vaccinator and air dry.
- Remove the air hose from the machine.
- Empty and wash ALL containers. Scrub all surfaces with a scrubbing pad and antibacterial detergent (dish washing liquid).
- Store all containers (exception – DI Water and Alcohol containers and alcohol) upside down to dry.
- Store distilled water and alcohol container on the device with ~0.5 gallon of alcohol (at least enough to cover the spigots), plug the lines back up and turn the spigots on.

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