DOCUMENT INFORMATION & PRINTING INSTRUCTIONS

ECN2981 – 6/19/2014 - JDC

NOTES:

1. THE FOLLOWING PAGES, NOT INCLUDING THIS COVER SHEET, TO BE PRINTED ON THE FRONT & BACK OF ONE PAGE, AND LAMINATED BASED ON THE FOLLOWING SPECIFICATIONS
2. ZOETIS PMS SPECIFIC 2-COLOR OR 4-COLOR PROCESS PRINTED ON 11” X 17” WITH SHORT-SIDE BLEEDS, 60# COVER-STOCK PAPER, DOUBLE-SIDED
3. LAMINATE INFORMATION: 5-MIL ENCAPSULATED WITH 1/4” BORDER, CORNERS ROUNDED
4. LABEL AND ARTWORK TO COMPLY WITH REQUIREMENTS OF ZOETIS LOGO STYLE GUIDE, 4-COLOR PMS MATCH OR PMS SPOT / GRADIENT COLORS

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF ZOETIS, INCORPORATED. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF ZOETIS, INCORPORATED IS PROHIBITED.

DAILY, WEEKLY, MONTHLY CLEANING PROCEDURES

DAILY CLEAN EMBREX® INOVOJECT® BIODEVICES AND EMBREX EGG REMOVER®

- **Disconnect** all power before washing device.
- Use an **antimicrobial** dishwashing detergent.
- Use scrub pads and brushes to remove organic matter.
- Use a long, narrow brush to clean tooling locator cups.
- DO NOT use high-pressure water or corrosive detergents.
- DO NOT spray water directly onto electrical components.
- DO NOT use abrasives on the ER or VS identifier/receiver lenses or HMI screens.
- DO NOT use abrasives or scrub pads on Interface.

DAILY ON-BOARD CLEANING STATION CONTAINERS

- Empty and rinse chemicals from all fluid containers.
- Spigots should be removed and disassembled for optimal cleaning
- Wash and scrub all containers inside and out including spigots.
- Rinse all soaps and detergents from containers and spigots.
- Store rinsed Sanitation and Cleaner containers upside down in a clean location to dry.
- DO NOT reuse old fluids.
- DO NOT leave fluids used for AM/PM cleaning in containers overnight.

NOTE: After the PM clean cycle has been completed and ALL containers are washed and rinsed, add 2 to 3 inches of alcohol to the DI water and Alcohol container and shake well. The containers should be hung on Inovoject with the lid secure and connected to the appropriate supply lines with spigot turned on. At start up, alcohol can be emptied from these two containers and replaced with appropriate fluids. Rinsing is not necessary.

DAILY SANITATION CONTAINERS

- Empty and rinse chemicals from containers.
- Scrub containers with soap and water; rinse out soap.
- Store containers upside down to dry in a clean location.
- DO NOT reuse old fluids.
- DO NOT leave fluids in containers overnight.
- DO NOT leave soap residues inside or outside of containers.

NOTE: Refill the 2-gallon alcohol container with fresh alcohol; store with pickup tube filter fully submerged in the container.

DAILY THAW BATH

- Empty and wash with soap and water.
- Dry and wipe down with alcohol.
- NOTE: Thaw bath can be covered or inverted for storage.
- DO NOT submerge electric thaw baths.
- DO NOT use high-pressure water or corrosive detergents.

DAILY ICE PACKS

- Clean daily with **antimicrobial** soap and rinse well.
- Dry and wipe down with alcohol.
- Return ice packs to the freezer.
- DO NOT use damaged or leaking ice packs.
- DO NOT leave ice packs out; return to freezer.

NOTE: Plastic/foil ice packs should be the ONLY type used. DO NOT use cloth ice packs.

DAILY VACCINE TRANSPORT COOLERS

- Empty and rinse, wash with **antimicrobial** soap and rinse.
- Dry and wipe down with alcohol using fresh paper towel.
- DO NOT store with fluids inside.
- DO NOT use for other purposes, e.g., food and beverages.

DAILY INTRATHERM BAGS

- Wash with soap and water; rinse well.
- Wash with other lab supplies.
- Hang to dry in clean location in the lab.
- Allow the bags to dry before reuse.
- DO NOT store bags on the device.
- DO NOT use the same brushes and pads used to wash the device.
- DO NOT use the bags for transportation of vaccine from the mixing room to the vaccinator.

WEEKLY ALCOHOL SPRAY BOTTLE AND NEEDLE CLEANING WIRE CONTAINER

- Empty and rinse the containers.
- Refill containers with fresh alcohol.
- DO NOT attempt to recycle used alcohol.
- DO NOT leave soap residues inside or outside of containers

WEEKLY LAB CLEANING

- Clean all fans (ceiling, free-standing, exhaust), including blades, motor, housing and louvers.
- Notify hatchery maintenance of defective air-handling units.
- Clean and disinfect tops of cabinets, counters, shelving and refrigerators.
- Clean and disinfect flooring.

MONTHLY LAB VENTILATION

- Clean lab air conditioning unit, change air filters for air conditioning unit, and clean air louvers inside and outside.
- Clean and disinfect refrigerator freezer (or stand-alone freezer).
- NOTE: Never use laboratory refrigerator/freezer for food or drink storage.

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