DOCUMENT INFORMATION & PRINTING INSTRUCTIONS

ECN2981 – 6/19/2014 - JDC

NOTES:

- 1. THE FOLLOWING PAGES, NOT INCLUDING THIS COVER SHEET, TO BE PRINTED ON THE FRONT & BACK OF ONE PAGE, AND LAMINATED BASED ON THE FOLLOWING SPECIFICATIONS
- 2. ZOETIS PMS SPECIFIC 2-COLOR OR 4-COLOR PROCESS PRINTED ON 11" X 17" WITH SHORT-SIDE BLEEDS, 60# COVER-STOCK PAPER, DOUBLE-SIDED
- 3. LAMINATE INFORMATION: 5-MIL ENCAPSULATED WITH 1/4" BORDER, CORNERS ROUNDED
- 4. LABEL AND ARTWORK TO COMPLY WITH REQUIREMENTS OF ZOETIS LOGO STYLE GUIDE, 4-COLOR PMS MATCH OR PMS SPOT / GRADIENT COLORS

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE
SOLE PROPERTY OF ZOETIS, INCORPORATED. ANY
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN
PERMISSION OF ZOETIS, INCORPORATED IS PROHIBITED.

HATCHERY CHECKLIST FOR EMBREX® INOVOJECT® BIODEVICES INSTALLATION AND USE

1. Adequate air supply for each area where Embrex® Biodevices will be used

- Ideally a ¾-inch ID supply line, rated at 300 PSI, running to the hatchery air compressor.
- A hatchery air compressor capable of 90 PSI at 15 CFM
- Clean, dry, filtered compressed air
- “Quick connects” supplied (should be installed within 50 feet of injection sites)

2. A clean, useable fresh water supply

- Ideally, hot water should be readily available in the hatchery

3. A storage area suitable for Embrex® BioDevices support equipment

- Well ventilated area that is not frequently fogged and away from chemical storage
- Low foot traffic

4. Dedicated electrical connections and supply cables readily available

- 30-amp designated GFI circuit breakers for electrical supply connections
- A fully functional ground connection in the hatchery electrical system
- 10-gauge 3 conductor (type SOW-A) electrical supply cables
- Supply connections should be installed within 50 feet of injection sites
- Strain relief for power connections (i.e., “twist-lock” style plugs)

5. Vaccine mixing supplies and preparation area available

- Vaccine
- Diluent and sterile dye
- Alcohol (70% isopropyl)
- Alcohol swabs and paper towels
- Syringes (20 ml)
- Needles (18 g x 1.5 inch)
- Distilled water

NOTE:
Ventilated vaccine prep area should be equipped with:

- Controlled and filtered airflow
- Thermostatically controlled vaccine thaw bath
- Adequate storage and work areas

6. Vaccine care and stabilization supplies, such as:

- Post-mixing transport cooler
- Crushed or cubed ice supply/cold packs (foil or plastic coated)
- Freezer

7. Cleaning and sanitizing chemicals

- Alcohol (70% isopropyl)
- Sterile saline or sterile water
- Non-foaming chlorinated cleaner
- Distilled water
- Inovotabs®
- Hand mixed sanitizer*
 - Chlorine bleach (5.25% concentration)
 - Sodium bromide
 - Citric acid
- Liquid dishwashing soap

* Directions for Hand mixed sanitizer can be found on PLACARD 3.

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