DOCUMENT INFORMATION & PRINTING INSTRUCTIONS

ECN2981 – 6/19/2014 - JDC

NOTES:

- 1. THE FOLLOWING PAGES, NOT INCLUDING THIS COVER SHEET, TO BE PRINTED ON THE FRONT & BACK OF ONE PAGE, AND LAMINATED BASED ON THE FOLLOWING SPECIFICATIONS
- 2. ZOETIS PMS SPECIFIC 2-COLOR OR 4-COLOR PROCESS PRINTED ON 11” X 17” WITH SHORT-SIDE BLEEDS, 60# COVER-STOCK PAPER, DOUBLE-SIDED
- 3. LAMINATE INFORMATION: 5-MIL ENCAPSULATED WITH 1/4” BORDER, CORNERS ROUNDED
- 4. LABEL AND ARTWORK TO COMPLY WITH REQUIREMENTS OF ZOETIS LOGO STYLE GUIDE, 4-COLOR PMS MATCH OR PMS SPOT / GRADIENT COLORS

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE
SOLE PROPERTY OF ZOETIS, INCORPORATED. ANY
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN
PERMISSION OF ZOETIS, INCORPORATED IS PROHIBITED.

EMBREX® INOVOJECT® BIODEVICES PROCESS CHECKLIST

This list contains the steps needed for verification and QC when using the Embrex® BioDevices. It is recommended that the vaccine mixer or crew leader perform this checklist DAILY at startup, after every vaccine bag change and throughout the day for early detection of operational problems.

INOVOJECT PROCESS

1. TOOLING “LOCATING TO EGG” PROPERLY ON EGGS
2. LOCATE AND REPLACE “BENT” NEEDLES
3. NO NEEDLE PUNCHING/NO BROKEN PUNCHES
4. NO “SLOW” INJECTION TOOLING CYLINDERS
5. NO “VISIBLE” VACCINE ON EGG SURFACES
6. NO FLUID LEAKS ON THE INJECTION PLATES
7. CLEAN / DRY IDENTIFIER “EYES”
8. USE QC PLATE PER 1% OF INJECTION CYCLES
9. NO “CLOGGED” VACUUM CUPS AT TRANSFER TABLE

1. Injection tooling is locating properly on top of the egg.

 - Tooling should align vertically with eggs so that the egg locator makes good contact with the egg surface.

NOTE: The most common cause for improper location is a misaligned pin stop.

 - Repositioning this pin stop will correct the tooling location. Contact a Zoetis Animal Health field service representative prior to making any adjustments.
2. No bent needles. If a needle is bent, change it.

 - Needles must be lowered for a visual inspection while the vaccine bag is replaced.

NOTE: Bent needles can be caused by misaligned tooling, improperly positioned pin stops, belt speeds or clamp tubing issues caused by fluids on injection head plates.
3. No needle punches. A visual inspection after injection is necessary in order to find needle punches.

 - The punch step is critical to ensure needle life, proper sanitation delivery and proper site of injection; a needle hole is visibly smaller than a normal punch hole.

NOTE: The 2 main causes of “needle punched” eggs are broken punches and fluid on the injection head and tooling (slipping tooling).
4. No slow injection tooling. Tooling speed should be checked visually after every vaccine bag change out.

 - To be effective, vaccine is dispensed only after a needle is fully extended, beneath the air cell membrane of the egg.

NOTE: Any needle not completely down before the vaccine pump dispenses is too slow.
5. No visible vaccine on the eggs. ANY visible vaccine on the egg surface should be investigated.

 - Common reasons for visible vaccine on the eggs are broken needles, partially clogged needles and tubes weeping.

NOTE: Weeping from the vaccine tubes can be a result of a loose vaccine pump cover or improper installation of tubes in the pump. For Vaccine Saver applications, pumps being mounted loose or duckbills sitting improperly under the pump may cause weeping.
6. Injection head is dry. Visually inspect injection head for fluids while waiting to hang a new bag.

 - Any fluid on the injection head can interfere with the process, contributing to needle punches and/or missed injections.
 - Any fluids leaking on the injection head will shorten the life of tooling and clamp tubing.
7. Identifiers are clean and dry.

 - Identifier assemblies particularly the lenses must remain clean and dry during operation
 - Always use a soft clean cloth

NOTE: NEVER use course paper towels or scrubbing pads to clean identifiers!
8. QC plate is run after every 1% of injection cycles.

 - If used properly, this test ensures that vaccine and sanitizer are being dispensed from every tooling.
 - All identifier assemblies should be cleaned using alcohol and paper towel during each QC cycle.**
9. Transfer head is operating properly.

 - Suction cups not releasing eggs into hatch basket will crush or drop eggs onto the next flat.
 - The delayed release of eggs is less obvious; cups should be observed DURING every transfer

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