## DOCUMENT INFORMATION & PRINTING INSTRUCTIONS

ECN2981 – 6/19/2014 - JDC

### NOTES:

1. THE FOLLOWING PAGES, NOT INCLUDING THIS COVER SHEET, TO BE PRINTED ON THE FRONT & BACK OF ONE PAGE, AND LAMINATED BASED ON THE FOLLOWING SPECIFICATIONS
2. ZOETIS PMS SPECIFIC 2-COLOR OR 4-COLOR PROCESS PRINTED ON 11" X 17" WITH SHORT-SIDE BLEEDS, 60# COVER-STOCK PAPER, DOUBLE-SIDED
3. LAMINATE INFORMATION: 5-MIL ENCAPSULATED WITH 1/4" BORDER, CORNERS ROUNDED
4. LABEL AND ARTWORK TO COMPLY WITH REQUIREMENTS OF ZOETIS LOGO STYLE GUIDE, 4-COLOR PMS MATCH OR PMS SPOT / GRADIENT COLORS

PROPRIETARY AND CONFIDENTIAL  
THE INFORMATION CONTAINED IN THIS DRAWING IS THE  
SOLE PROPERTY OF ZOETIS, INCORPORATED. ANY  
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN  
PERMISSION OF ZOETIS, INCORPORATED IS PROHIBITED.

## INOVOTABS® SANITIZER PREPARATION

1. Prepare a fresh solution at the **START** of each day.
2. Ensure that the container is washed, rinsed and free of any sediment or organic matter before use.
3. Use **ONLY** potable water (water suitable for drinking) to mix Inovotabs® solution.
4. Add 2 Inovotabs tablets per gallon (1 tablet/2 quarts or 2 liters) of water.
  - Warm water may be used (75-104 F/23-40 C).

**CAUTION: The solution must not be hot when the machine is injecting eggs.**

5. Wait 15 minutes for Inovotabs to dissolve; gently agitate before use.
  - Ensure that all Inovotabs have been completely dissolved prior to their use in Embrex® BioDevices.
6. Cover/Close the containers with a cap after mixing.
7. Connect containers to Inovoject® and follow normal procedures for operation.

During the flushing of the sanitation system, it is imperative that water be flushed between sanitizer and alcohol for a minimum of 2 gallons; this can be done during the AM/PM vaccine system cleaning cycle.

## INOVOTABS SAFETY INFORMATION

Follow all appropriate safety guidelines when handling tablets.

Please read all precaution information on label located on the side of the product container.

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